# Menuett®

## 013697

### **SV FEBERTERMOMETER**

#### BRUKSANVISNING

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## SAFETY INSTRUCTIONS

Read all the instructions and safety instructions carefully before use and save them for future reference.

### WARNING!

- **Suffocation risk: the thermometer casing and battery can fasten in the throat and lead to suffocation and death. Do not allow children to use the product without supervision.**
- **The product is not intended for body temperature measuring in the auditory canal. The product is only intended for body temperature measuring in the rectum, oral cavity and armpit.**
- **Do not expose the product battery to naked flames or high temperatures – risk of explosion.**
- **Remove the battery if the product is not going to be used for some time.**
- **The product cannot replace consultation with a doctor. Consult a doctor if you need to interpret the readings or are concerned about your health status. It can be dangerous to diagnose yourself and could cause a deterioration in your health.**
- **Do not use the product if it is wet – risk of unreliable result.**
- **Do not bite the product – risk of personal injury and damage to the product.**
- **Do not dismantle the product or attempt to modify, or repair it. Tampering with the product can make it unreliable. Disinfect the product after use, especially if it has been used by several persons. For use in the rectum, insert the product carefully. If experiencing pain, stop the insertion – risk of personal injury.**
- **Do not use the product orally after use in the rectum.**
- **Do not use the product orally for children less than two years old.**
- **If the product has been stored at a temperature outside the interval of 5 to 40°C, allow the temperature of the product to equalise for at least 15 minutes in a room where the temperature is within the interval 5 to 40°C before use.**

### NOTE:

- **If the product is used with a probe cover, the temperature reading can deviate by 0.1°C (0.2°F) from the correct value.**
- **Product performance can be reduced by:**
  - Using at temperatures or humidity outside the interval specified by the manufacturer.
  - Storage at temperatures or humidity outside the interval specified by the manufacturer.
  - Knocks and blows, or other mechanical stress.
  - Patient temperature lower than the ambient temperature.
- **Portable and mobile radio communication equipment can affect the product. Special conditions concerning electromagnetic interference must be complied with, see section “Electromagnetic compatibility”.**

## SYMBOLS

Read the instructions.

	Approved in accordance with the relevant directives.
	Recycle discarded product in accordance with local regulations.
	European representative.
	Direct voltage.
	Batch code.
	Manufacturer.
	Transport and storage temperature: $-20$ to $55^{\circ}\text{C}$ ( $-4$ to $131^{\circ}\text{F}$ ).
	BF type, patent applied for part.
	Medical Device.

## TECHNICAL DATA

Battery type	Button cell battery
Battery	1.5 V (1 x LR41, SR41 or UCC 392)
Battery life	Approx. 200 h continuously or 3 measurements per day for 1 year
Size	L123 x W18 x H9 mm
Weight	approx. 10 g
Expected life span	3 years
Display	LCD, 3 1/2 digits
Duty cycle	Direct display
Ambient temperature	5 to $40^{\circ}\text{C}$ , relative humidity 15 to 95% RH Atmospheric pressure 700 to 1060 hPa
Storage temperature	$-20$ to $55^{\circ}\text{C}$ , relative humidity 15 to 95% RH Atmospheric pressure 700 to 1060 hPa
Protection rating	IP22
Classification	BF type
Measurement range	$32.0$ to $42.9^{\circ}\text{C}$ ( $90.0$ to $109.9^{\circ}\text{F}$ ) (temperature unit $^{\circ}\text{C}$ or $^{\circ}\text{F}$ is preselected on delivery)

Accuracy  $\pm 0.1\text{ C}$  ( $0.2^{\circ}\text{F}$ ) in measurement value interval  $35.5$  to  $42.0^{\circ}\text{C}$  ( $95.9$  to  $107.6^{\circ}\text{F}$ ) at ambient temperature  $18$  to  $28^{\circ}\text{C}$ .  $\pm 0.2^{\circ}\text{C}$  ( $\pm 0.4^{\circ}\text{F}$ ) outside these temperature intervals.

The product complies with requirements in MDR(EU)2017/745, approved by notified body no. 0123.

### Manufacturer

JOYTECH HEALTHCARE CO. LTD. No. 365, Wuzhou Road, Yuhang Economic Development Zone, Hangzhou city, 311100 Zhejiang, China.

### European representative

Shanghai International Holding Corp. GmbH (Europe)  
Eiffestrasse 80, 20537 Hamburg, Germany

## DESCRIPTION

The product is a digital fever thermometer for quick and accurate measuring of body temperature. The product is a fever thermometer for oral, rectal and auxiliary measuring of body temperature, and is intended to be used multiple times in clinics and in private homes. The fever thermometer can be used for people of all ages. For best results and product life span, read these instructions carefully before use and save them for future reference.

The product complies with the following standards:

- SS-EN 12470-3 Medical thermometers - Part 3: Electrical thermometers (direct display and calculation) with maximum function.
- ISO 80601-2-56 Electrical equipment for medical use - Part 2-56: Specific requirements for basic safety and essential performance for medical thermometers for measuring body temperature.
- SS-EN 60601-1-11 Electrical equipment for medical use - Part 1-11: General requirements concerning safety and essential performance – Supplementary standard for equipment and systems for use in domestic care environments. The product also complies with EN 60601-1-2 (EMC) and IEC/EN60601-1 (Safety).
- The manufacturer is certified in accordance with ISO 13485.

## PACKAGE CONTENTS

- Thermometer
- Instructions
- Storage case

## PARTS

1. Probe
2. Display
3. Power switch
4. Battery cover

FIG. 1

## USE

1. Press the power switch button beside the display. An audio prompt is heard and the display first shows **188.8<sup>°E</sup>** and thereafter the last registered temperature. After showing the self-test temperature the fever thermometer is ready to use.
2. Place the thermometer at the required measuring point – oral, rectal or auxiliary.
  - Oral measuring (oral cavity): Put the thermometer under the tongue, in one of positions marked “V” in figure 2. Close the mouth and breath evenly and calmly through the nose, so that the air flow when breathing in and out does not affect the measurement. The oral temperature of a healthy person is normally between 35.7 and 37.3°C (96.3 and 99.1°F).

FIG. 2

- Rectal measuring (rectum): Apply vaseline to the tip of the probe to simplify insertion. Carefully insert the probe about 1 cm into the rectum. The rectal temperature of a healthy person is normally between 36.2 and 37.7°C (97.2 and 99.9°F).
  - Auxiliary measuring (armpit): Wipe the armpit dry. Put the probe in the armpit and hold the arm firmly pressed against the side.  
This method of measuring gives a result that is not medically reliable. One of the other methods should be used for accurate and medically reliable results. The auxiliary temperature of a healthy person is normally between 35.2 and 36.7°C (95.4 and 98.1°F).
3. The degree symbol flashes during the measuring. The degree symbol stops flashing when the measurement is complete and an intermittent audio prompt is given for about 10 seconds. The measured temperature is shown on the display at the same time.  
The measurement must not be interrupted before audio prompt indicates that the measurement is done. The thermometer continues to measure the temperature even after the audio prompt has indicated that the measurement is done. This can be used to improve the accuracy and reliability of the measurement result: continue holding the probe at the measuring point for at least 30 seconds, preferably longer – about two minutes for oral or rectal measuring, about five minutes for auxiliary measuring.
  4. The thermometer switches off automatically after about 10 minutes of inactivity. To save the batteries the thermometer can also be switched off manually with the power switch as soon as the measurement is ready.

### NOTE:

**The audio prompt consists of a uniform series of beeps “tut --- tut --- tut --- tut”. When the measured temperature is 37.8°C (100°F) or higher the beeps come more quickly and groupwise “tut - tut - tut ----- tut - tut - tut ----- tut - tut - tut”.**

## MAINTENANCE

### CLEANING

- Wipe the surface of the thermometer with a soft, clean cloth.
- Stubborn stains can be removed with a soft cloth moistened with a mild, non abrasive detergent. Wipe dry with a soft, dry cloth. Disinfect with a 75% solution of ethanol or isopropanol.

**WARNING!**

- **Do not use detergent that contains benzine, thinner, petrol or other strong solvents.**
- **Do not attempt to disinfect the probe (measuring tip) by dipping it in alcohol or hot water (50°C or hotter).**
- **Do not wash with ultrasound.**

**REPLACING THE BATTERY**

1. The battery should be replaced when the symbol  is shown in the bottom right corner.
2. Remove the battery cover by pushing it back.

**FIG. 3**

3. Carefully pull out the circuit board on which the battery sits about one cm.

**FIG. 4**

4. Remove the used battery with a pen or similar sharp implement. Insert a new 1.5 V button cell battery LR41, SR41, UCC392 or the equivalent. Insert the battery with the polarity marking “+” facing up.

**FIG. 5**

5. Push in the circuit board and replace the battery cover.

**CALIBRATION**

The fever thermometer is supplied calibrated from the factory. No regular recalibration is necessary when used in accordance with these instructions. It is recommended to calibrate the fever thermometer every two years, or if its clinical accuracy is doubtful. Start the thermometer and dip the probe in a water bath and check if the temperature reading differs unacceptably from a laboratory approved reference thermometer.

The thermometer should then be submitted to the manufacturer or retailer for recalibration.

The above recommendation does not replace any legal requirements. Legislation, rules and regulations, and guidelines for the use, functionality and accuracy of the fever thermometer should always be complied with.

**TROUBLESHOOTING**

Error message	Problem	Action
	Measured temperature is less than 32.0°C (90.0°F).	Switch off the thermometer. Wait one minute, start the thermometer and attempt another measurement with good surface contact.
	Measured temperature exceeds 42.9°C (109.9°F).	Switch off the thermometer. Wait one minute, start the thermometer and attempt another measurement with good surface contact.

	The product is not working correctly.	Remove the battery, wait one minute, replace the battery, and try to use the thermometer again. Contact the retailer if the fault persists.
	Low battery charge.	Replace battery.

## ELECTROMAGNETIC COMPATIBILITY

The product complies with the EMC requirements in accordance with the international standard IEC 60601-1-2 if the specified conditions in the table below are met. The product is an electrical medical product and as such is subject to special requirements concerning electromagnetic compatibility, which requirements must be indicated in the instructions for the product. Portable and mobile radio communication equipment can affect the product. Do not use the product together with unapproved accessories – this can have a negative effect on the product and change its electromagnetic compatibility. Do not use the product in the immediate vicinity of electrical equipment or between different electrical equipment.

## GUIDELINES & MANUFACTURER'S CERTIFICATE CONCERNING ELECTROMAGNETIC RADIATION

The product is intended to be used in electromagnetic environments in accordance with the specifications below. It is incumbent on the customer or the person using the product to ensure that the electromagnetic environment in which the product is used complies with these specifications.

Emission test	Compliance	Guidelines for electromagnetic environment
Radio frequency radiation CISPR 11	Group 1	Radio frequency energy is only used in the product for the internal functions in the product. The emitted radio frequency radiation is therefore very low and the probability of it causing interference in nearby electronic equipment is low.
Radio frequency radiation CISPR 11	Class B	Class B equipment can be used for domestic purposes and in facilities directly connected to low voltage networks in residential buildings.
Overtones IEC 61000-3-2	Not applicable	
Voltage fluctuation/flicker IEC 61000-3-3	Not applicable	

## GUIDELINES & MANUFACTURER'S CERTIFICATE CONCERNING ELECTROMAGNETIC IMMUNITY

The product is intended to be used in electromagnetic environments in accordance with the specifications below. It is incumbent on the customer or the person using the product to ensure that the electromagnetic environment in which the product is used complies with these specifications.

Immunity test	Test level as per IEC 60601	Compliance level	Guidelines for electromagnetic environment
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 8$ kV line carrier $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV airborne	$\pm 8$ kV line carrier $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV airborne	The floor surface should consist of wood, concrete or clinker. In rooms with floors covered with synthetic material the relative air humidity should be at least 30%.
Electrostatic rapid transient/burst IEC 61000-4-4	$\pm 2$ kV for power lines $\pm 1$ kV for incoming and outgoing lines	Not applicable	
Overvoltage pulse IEC 61000-4-5	$\pm 1$ kV differential mode $\pm 2$ kV common mode	Not applicable	
Undervoltage pulses, short-term interruption and voltage variations on input lines. IEC 61000-4-11	< 5% OUT (> 95% reduction of OUT) during 0.5 period 40% OUT (60% reduction of OUT) during 5 periods 70% OUT (30% reduction of OUT) during 25 periods < 5% OUT (> 95% reduction of OUT) during 5 seconds	Not applicable	
Net frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50/60 Hz	30 A/m, 50/60 Hz	The magnetic field of the net frequency should be at the levels that characterise a normal place of use in a typical commercial environment or hospital environment.

## GUIDELINES & MANUFACTURER'S CERTIFICATE CONCERNING ELECTROMAGNETIC IMMUNITY

The product is intended to be used in electromagnetic environments in accordance with the specifications below.

It is incumbent on the customer or the person using the product to ensure that the electromagnetic environment in which the product is used complies with these specifications.

Immunity test	Test level as per IEC 60601	Compliance level	Guidelines for electromagnetic environment
Line carried radio frequency interference IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not applicable	Portable and mobile radio communication equipment must not be used nearer to any part of the product or its leads than the recommended minimum distance in accordance with the equation applicable for the transmitter in question.
Emitted radio frequency interference IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m	Recommended minimum distance $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz
Radio communication equipment IEC 61000-4-3	380 MHz, 27 V/m 450 MHz, 28 V/m 710 MHz, 745 MHz, 780 MHz 9 V/m 810 MHz, 870 MHz, 930 MHz 28 V/m 1720 MHz, 1845 MHz, 1970 MHz 28 V/m 2450 MHz, 28 V/m 5240 MHz, 5500 MHz, 5785 MHz 9 V/m	380 MHz, 27 V/m 450 MHz, 28 V/m 710 MHz, 745 MHz, 780 MHz 9 V/m 810 MHz, 870 MHz, 930 MHz 28 V/m 1720 MHz, 1845 MHz, 1970 MHz 28 V/m 2450 MHz, 28 V/m 5240 MHz, 5500 MHz, 5785 MHz 9 V/m	$d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.7 GHz where P is the maximum output specified by the transmitter manufacturer in watts (W) and d is the recommended minimum distance in metres (m). The field strength from a permanent radio transmitter, determined by electromagnetic measurement in situ, should be less than the specified compliance level for the respective frequency range. Interference can occur in the vicinity of equipment marked with the symbol: 

## RECOMMENDED MINIMUM DISTANCE BETWEEN THE PRODUCT & PORTABLE & MOBILE RADIO COMMUNICATION EQUIPMENT

The product is intended to be used in electromagnetic environments in which radiated radio frequency interference is controlled and limited. The customer or person using the product can reduce the risk of electromagnetic interference by maintaining the recommended minimum distance, given below on the basis of radio transmitted power, between the product and portable and mobile radio communication equipment.

## RECOMMENDED MINIMUM DISTANCE BETWEEN THE PRODUCT & PORTABLE & MOBILE RADIO COMMUNICATION EQUIPMENT

Rated output power P in watts (W) of radio transmitter	Minimum distance d between product and radio transmitter in metres (m)	
	80 MHz till 800 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz till 2,7 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0,01	0,12	0,23
0,1	0,38	0,73
1	1,2	2,3
10	3,8	7,3
100	12	23

For radio transmitters where the rated power P is not given in the table values the minimum distance d can be calculated with the applicable equation for the transmitter frequency, whereby the transmitting power P is given in W.

### NOTE:

- The minimum distance for the higher frequency range is applicable at 80 MHz and 800 MHz.**
- These guidelines do not apply in all situations. Electromagnetic radiation can be affected by absorption in and reflection against building structures and other objects, such as people.**

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

### Hersteller

JOYTECH HEALTHCARE CO.LTD.No. 365, Wuzhou Road, Yuhang Economic Development Zone, Hangzhou city, 311100 Zhejiang, China.

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Shanghai International Holding Corp. GmbH (Europa) Eiffestrasse 80, 20537 Hamburg, Deutschland

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© Jula AB