EN FOOT WARMERS
OPERATING INSTRUCTIONS
Original instructions

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Jula AB reserves the right to make changes to the product. Jula AB claims copyright on this documentation. It is not allowed to modify or alter this documentation in any way and the manual shall be printed and used as it is in relation to the product. For the latest version of operating instructions, refer to the Jula website.

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
1 Introduction

1.1 The product

The product is a single use warmer for feet that gets activated by contact with air. The product keeps a temperature of 40°–65°C for a minimum of 5 hours.

1.2 Dedicated use

The product can be used to keep feet and toes warm, or to give temporary relief from small aches and pains of muscles and joints.

1.3 Symbols

	Read the operating instructions carefully and make sure that you understand the instructions before you use the product. Save the instructions for future reference.
	This product complies with applicable EU directives and regulations.
	Discarded product must be recycled according to current regulations.
	Caution!
	Keep dry.
	Keep away from sunlight.
	Do not use the product if damaged, refer to instructions for use.
	Do not reuse the product.
	The product is a medical device.

2 Safety

2.1 Safety definitions

WARNING! If you do not obey these instructions, there is a risk of death or injury.

CAUTION! If you do not obey these instructions, there is a risk of damage to the product, other materials or the adjacent area.

NOTE! Information that is necessary in a given situation.

2.2 Safety instructions for operation

NOTE! Report all serious incidents related to the device to the manufacturer and to the competent authority in the Member State where the user and, or patient is established.

• Read and obey the warning instructions before operation.

• Do not let children or persons who do not know the product use it.

• Do not use the product if it is damaged.

• Do not make changes to the product.

• Use the product only for its specified function.

• Do not use the product directly on your skin.

• Do not use the product on children, or persons who can not remove the product themselves.

• Do not use the product when you sleep.

• Do not use on areas of skin that are sensitive to heat.

• Refer to a physician or doctor before use if you have circulatory problems, rheumatoid arthritis, nerve damage, sensory impairment, skin conditions, are diabetic or pregnant.

• Do not put the product in your mouth or eyes. The product is for external use only.

3 Operation

NOTE! The product activates when it get in contact with air. Keep the product in its packaging until you use it.

NOTE! If the product become too hot, put it in cloth.

NOTE! Use the product only one time. Discard the product after use.

3.1 To operate the product

1 Open the outer packaging.

2 Shake the product for 30 seconds.

3 Remove the protective cover from the product.

4 Attach the adhesive side of the product to the outer surface of your sock at the the specified area of your foot.

4 Discard

• Make sure that you follow local regulations when you discard the product. Do not burn the product.

5 Technical data

Specification	Value
Size	9 cm x 7 cm
Heat duration	5 hours
Temperature	40–65°C
Ingredients	Iron powder, vermiculite, water, activated carbon, salt

RFX+CARE Manufacturing Co., Ltd.
No.599 Wuxie Road, Yuecheng District, Shaoxing,
Zhejiang, China 312000

RFX+Care International A/S
Gammel Strandvej 22B 2990 Nivå, Denmark
www.rfxcare.com

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Jula AB, Box 363,532 24 SKARA, Sweden-www.jula.com